

UNIVERSITY OF ILLINOIS
AT URBANA-CHAMPAIGN

Physics 403. Modern Physics Laboratory

Fall 2020

Eugene V. Colla, Virginia O. Lorenz
COVID-19 - hybrid in-person-online version



illinois.edu

Physics 403 Modern Physics Laboratory

Fall 2020 Teaching Team



Eugene




Gina



April

April Townsend
aprilmf2@illinois.edu



Andrew

Andrew Calhoun
ajc7@illinois.edu



Abid

Abid Khan
aakhan3@illinois.edu



Albur


Hassan, Albur
ahassan4@illinois.edu

Instructors:

Eugene V Colla
kolla@illinos.edu

Virginia (Gina) Lorenz
vlorenz@illinois.edu

Support from Paul Kwiat Team




Jack

**Laboratory
Specialist:**
Jack Boparai
jboparai@illinois.edu



Andrew

Andrew Conrad
aconrad5@illinois.edu



Sam

Samantha Isaac
isaac5@illinois.edu



Spencer

Spencer Johnson
sjj3@illinois.edu

Outline

- I. Goals of the course**
- II. Teamwork / grades / expectations from you**
- III. Syllabus and schedule**
- IV. Your working mode**
 - In class and “after hours” access
 - Safety, Responsibility
 - Home and away computing
- V. Take a Lab tour (only video)!**
- VI. Let’s get started**
 - electronic logbooks



Course Goals. Primary goals:

- **Learn how to “do” research**

- ✓ **Each project is a mini-research project**

- ✓ **How are experiments actually carried out ?**

The procedures aren't all written out

The questions are not in the back of the chapter

The answers are not in the back of the book

You will have to learn to guide your own activities

- ✓ **Use of modern tools and modern analysis and data-recording techniques**



Course Goals. Primary goals:

- **Learn how to document your work**
 - **Online - electronic logbook ***
 - **Online – saving data and projects in student area on server**
 - **Using traditional paper logbooks**
 - **Making an analysis report**
 - **Writing formal reports**
 - **Presenting your findings orally (online)**



Course Goals. Secondary goals:

- **Learn some modern physics**
 - **Many experiments were once awarded by Nobel-prize**
 - **They touch on important themes in the development of modern physics**
 - **Some will provide additional insight to understand advanced courses you have taken**



The Experiments. Three main groups





- **Nuclear / Particle (NP)**
- **Atomic / Molecular / Optics (AMO)**
- **Condensed Matter (CM)**

You will do the experiment from all these groups



The Experiments

- **Nuclear / Particle (NP)**
 - Alpha particle range in gasses
 - γ - γ correlation experiment
 - γ – spectroscopy
 - Mössbauer spectroscopy



The Experiments

- **Nuclear / Particle (NP)**



- Cosmic ray muons:

Lifetime, capture rate, magnetic moment

- Angular distribution of cosmic rays

- γ – spectroscopy




- Mössbauer spectroscopy (new)



The Experiments

Atomic/Molecular/Optics (AMO)




- Berry's phase
- Quantum erasure
- Quantum Entanglement



The Experiments

Atomic/Molecular/Optics (AMO)




- Optical pumping of rubidium gas
- Fluorescence spectroscopy



The Experiments





- Condensed Matter (CM)
 - Superconductivity
 - Tunneling in superconductors
 - 2nd sound in ⁴He superfluid

state




The Experiments

- **Condensed Matter (CM)**
 - Ferroelectrics and ferroelectric phase transition
 - Pulsed NMR
 - Calibration of temperature sensors





The Experiments

- **Condensed Matter (CM)**
 - Special Tools:
 - Vacuum film deposition
 - Atomic Force Microscope
 - Polarizing microscope



The “manuals”

- Many are just guides
- A only few purchased experiments have “real” manuals
- We serve as your guides ... like real research ... yes, we will do it in “online” mode too. We have prepared the materials explaining how to do the experiments and data analysis and you can find all these materials and examples of data analysis in folder in common drive.



TEACH
SPIN
Instruments Designed for Teaching

OPTICAL PUMPING
OF RUBIDIUM
OP1-A



Outline

I. Goals of the course

II. Teamwork / grades / expectations from you

III. Syllabus and schedule

IV. Your working mode

In class and “after hours” access

Safety, Responsibility

Home and away computing

V. Take a Lab tour !

VI. Let's get started

electronic logbooks

digital scopes



Grading: Distribution of “740” points

ASSIGNMENT	Points
Expt. documentation: elog reports, shift summaries, plot quality; paper logbooks	120 Total 60 / cycle
Formal reports: physics case, quality of results, depth of analysis, conclusions	400 Total 100 / report
1st Oral report: motivation, organization of presentation; fielding questions	100 Total
Final Oral Presentation \equiv Final Exam	120
Total Effective point total will be	740


The grading scale will be a percentage out of “740” :

Letter grading scale is approximately **97% = A+**, **93% = A**, **90% = A-**, **87% = B+**, **83% = B**, **80% = B-**, etc

You can **RESUBMIT one lab report** to improve your grade
(deadline for resubmissions and for report #4 **Dec. 8th 2020**)



Grading: a piece of history and analysis of the results



Submission of Lab-Reports

- Due dates as on syllabus at midnight
- The reports should be uploaded to the server:
- <https://my.physics.illinois.edu/courses/upload/>
- Accepted MS-Word or PDF
- For orals – MS-PowerPoint* or PDF

** preferable*



Absences

- If you are sick, **let Eugene know by email (kolla@Illinois.edu)**. Don't come in and get others sick. We are working side-by-side in a close environment for many hours.
- You can “make up” the time with arrangements and you can have access to the rooms. We will be accommodating.



Absences. Excuse Policy.

- You can be excused from **only one** missed assignment, and only if you **provide medical or any other acceptable documentation**¹.
- If the excused you have **missed the oral presentation** (oral #1), you have to discuss this with us and we will arrange the date for your oral talk.
- The Final Oral **cannot be excused**, as it is equivalent to a final exam. You **cannot pass the course without credit for this assignment**²

1. Student Code: <https://studentcode.illinois.edu/article1/part5/1-501/>

2. Ibid: <https://studentcode.illinois.edu/article3/part2/3-201/>



Late Reports

- **Policy for late reports**

- You can have **ONE “late ticket”** for a **“free”** delay of up to **3 business days**, but you must tell us you are using the ticket
- Reports are due at midnight on the date shown on the syllabus. After that we will charge:
 - 5 points for up to 1 week late. 10 points for up to 2 weeks late.
 - After that, it's too late.



C1-Ex1(2.07.18)

Outline

- I. Goals of the course
- II. Teamwork / grades / expectations from you
- III. Syllabus and schedule**
- IV. Your working mode
 - In class and “after hours” access
 - Safety, Responsibility
 - Home and away computing
- V. Take a Lab tour !
- VI. Let’s get started
 - electronic logbooks
 - digital scopes



Syllabus

Cycles

	Date	Day	Activity	Comment	Note
1	8/25	Tues	Orientation	About Phy403	online
2	8/27	Thurs	Cycle 1-1		
3	9/01	Tues	Cycle 1-2	OriginPro Intro/Root	
4	9/03	Thurs	Cycle 1-3	Elog Comments	
5	9/08	Tues	Cycle 1-4	Written Reports	
6	9/10	Thurs	Cycle 1-5		
7	9/15	Tues	Cycle 1-6	Error analysis	
8	9/17	Thurs	Cycle 1-7		
9	9/22	Tues	Cycle 1-8	Oral Reports/Talks	C1-Ex1 (9.23.20)
10	9/24	Thurs	Cycle 1-9		
11	9/29	Tues	Cycle 1-10	Optical spectroscopy	
12	10/01	Thurs	Cycle 1-11		
13	10/06	Tues	Cycle 1-12	Ferroelectricity	
14	10/08	Thurs	Cycle 2-1		Rotate
15	10/13	Tues		ORALS 1	online
16	10/15	Thurs			online
17	10/20	Tues	Cycle 2-2	High Energy Physics	C1-Ex2 (10.21.20)
18	10/22	Thurs	Cycle 2-3		
19	10/27	Tues	Cycle 2-4	Noise (mw)	
20	10/29	Thurs	Cycle 2-5		
	11/03			Election Day	
21	11/05	Thurs	Cycle 2-6	Lock-in Amps and FT	
22	11/10	Tues	Cycle 2-7		
23	11/12	Thurs	Cycle 2-8	Entanglement	C2-Ex1 (11.13.20)
24	11/17	Tues	Cycle 2-9		
25	11/19	Thurs	Cycle 2-10	Measuring Temp	
				Thanksgiving Break	
26	12/01	Tues	Cycle 3-11		online
27	12/03	Thurs	Cycle 3-12	To be announced	online
28	12/08	Tues		Final Orals #1	online
	12/10			READING DAY	C2-Ex2 (12.10.20)
29	12/15	Thurs		Final Orals #1	online

* Lecture topics are subject to change



	NP	CM	Atomic + CM	Optics
	A. Cosmic Muon Stand i. Muon lifetime ii. Capture rate iii. Magnetic moment B. Alpha range C. Gamma Gamma D. Muon telescope E. Mössbauer spectroscopy	A. Ferro 1 B. Ferro 2 (imaging) C. 2 nd sound of ⁴ He D. Hysteresis loops E. Tunneling F. T calibration	A. Optical pumping B. Superconductivity C. Mutual inductance D. pNMR	A. Quantum Table i. Berry's phase ii. Quantum erasure iii. Entanglement B. Fluorescence spectroscopy C. AFM
	Virginia, April	Eugene	Eugene, Albur, Andrew	Abid, TAs from Kwiat Lab
C1-1	1-2, 3-4, 5-6	7-8, 9-10, 11-12	13-14, 15-16, 17-18	19-20, 21-22, 23-24
C1-2	19-20, 21-22, 23-24	1-2, 3-4, 5-6	7-8, 9-10, 11-12	13-14, 15-16, 17-18
C2-1	14-16, 15-17, 18-13	20-22, 21-23, 24-19	2-4, 3-5, 6-1	8-10, 9-11, 12-7
C2-2	8-10, 9-11, 12-7	14-16, 15-17, 18-13	20-22, 21-23, 24-19	2-4, 3-5, 6-1



Cycle	#	Experiment
C1-1	1-2	Cosmic ray muons
	3-4	Alpha range
	5-6	Gamma-gamma
	7-8	Ferro 1
	9-10	Ferro 3
	11-12	Second Sound
	13-14	NMR
	15-16	Superconductivity
	17-18	Optical Pumping
	19-20	Fluorescence
	21-22	Quantum Optics
	23-24	Quantum Optics
C1-2	19-20	Cosmic ray muons
	21-22	Alpha range
	23-24	Mössbauer spectroscopy
	1-2	Ferro 1
	3-4	Tunneling
	5-6	Ferro 2
	7-8	NMR
	9-10	Superconductivity
	11-12	Optical Pumping
	13-14	Fluorescence
	15-16	AFM
	17-18	Quantum Optics



Assignment of experiments

2 cycles with 2 experiments

→ teams change after cycle

→ joint team reports and elogs but oral


presentations will be done by each

student personally



Spring 2019 Orals Physics 403

After 2 experiments (1 cycle) we will have oral session. The topic of the presentation will be chosen from the experiments done in this cycle.



Outline

- I. Goals of the course
- II. Teamwork / grades / expectations from you
- III. Syllabus and schedule
- IV. Your working mode**
 - In class and “after hours” access
 - Safety, Responsibility
 - Home and away computing
- V. Take a Lab tour !
- VI. Let’s get started
 - electronic logbooks
 - digital scopes



Fall 2020 working mode*

This semester we plan to work in hybrid mode:

Each experiment will take **six lab days** and you are going to work in **team of two**.

To keep in the lab **only 12 P403 students** we will share the granting the access to the lab only half of the class at the same time.

It means that one lab day the **partner no1** will in Lab and **partner no2** – online.

Next lab day you will be swapped – **no2** in lab and **no1** online.

Finally each student will be able to work in lab during 3 days (50% of total time).



** Subject to change dependable on COVID19 situation and “recommendations” provided by Uofl administration.*



Fall 2020 working mode.

Each team works together and will have the common grades for report and elogs. It is possible some variations of the prepared team schedule dependable on presence of teammates. The P403 lab **works in real time** according the course schedule and both partners should work on experiments during the **whole lab time in person or online.**



Lab Access



Use Your ID Card to Access the Lab

You can access the Lab not only on “Lab days”

Late time rules:

~~You can stay in the Lab until 8pm but need to
work with partner~~

~~After 8pm and on weekend days – you have to discuss
this schedule with your instructor and in general it is
preferable to avoid working after 8 pm and on week~~

These will be not the options for Fall 2020 semester



Safety is your responsibility !

Hazards: *high voltage, radioactive sources,
cryogens, chemical materials, high pressure*

In class work and “after-hours” access work requires
responsible conduct with regards to

(I) safety/hazards and with

(II) equipment

Discuss potential hazards at the beginning of each
experiment with an instructor or TA

When in doubt stop and ask



Follow Directly the Recommendations of Safety Working

<https://www.drs.illinois.edu/>

RESEARCH SAFETY

Accident Response ▾ DRS Safety Programs ▾ Training ▾ Waste Management ▾ Safety Library ▾

My Campus User Login

(Material) Safety Data Sheets

Find safety data sheets for material you work with in your lab.

NEWS AND ANNOUNCEMENTS [VIEW ARCHIVE »](#)

Laser Registration and Management
9/23/2018
The Division of Research Safety has added a tool to their website to allow laser users to manage their laser registrations and inventory on-line.

New Tier 1 Select Agent
9/23/2018
As of 9/14/16, the CDC/HHS has added Bacillus cereus Biovar anthracis as a Tier 1 select agent under 42 CFR Part 73.

Laser Safety Eyewear Warning
7/6/2018
Filters not matching specifications on packaging

RESPONSIBILITIES

I work in a lab

I supervise a lab

My work takes me into labs

Follow Directly the Recommendations of Safety Working

Related Units @ Illinois Questions?

Division of
RESEARCH SAFETY

Accident Response ▾ DRS Safety Programs ▾ Training ▾ Waste Management ▾ Safety Library ▾

Chemical Waste Collection and Storage

Before generating chemical waste, the researcher should determine how it will be collected and stored and obtain the necessary equipment (containers, labels) in advance. The choice of procedures depends on the type of waste and its final disposition. This section explains how to determine the final disposition of waste, select the appropriate waste container, and store waste in the lab or work area. It also suggests waste minimization strategies.

Determining How to Dispose of a Chemical Waste

The final disposition of a chemical waste is determined by the answers to a series of questions:

Step 1. Is the waste [Contaminated Debris](#) (glassware, paper towels, clean-up materials), or is it a chemical or chemical mixture?
If it is contaminated debris: Go to Step 5.
If it is a chemical or chemical mixture: Go to Step 2.

Step 2. Is the chemical a DEA (Drug Enforcement Agency) controlled substance? (Refer to the [DEA list controlled substances](#).)
Yes: Refer to the [DEA Controlled Substances Guide](#) for disposal procedures.
No: Go to Step 3.

Step 3. Is the chemical a solid (not liquid or gas)?
Yes: Collect and store the waste as described in the waste container and storage guidelines listed below and dispose of it through the Division of Research Safety (DRS) chemical waste disposal program. See the section [Procedures for Requesting Chemical Waste Disposal](#) for the disposal procedures. (No solid chemical waste, hazardous or non-hazardous, should be placed in the regular trash.)
No: Go to Step 4.

Step 4. Is the chemical a liquid non-hazardous waste as listed in the section [Liquid Non-Hazardous Chemical Waste Disposal](#)?
Yes: The chemical may be poured down the sanitary sewer (sink drain) with copious amounts of water.
No: Collect and store the waste as described in the waste container and storage guidelines listed below, and dispose of it through the DRS chemical waste disposal program. See the section [Procedures for Requesting Chemical Waste Disposal](#) for the disposal procedures.

Step 5. Is the contaminated debris laboratory glassware (broken and unbroken)?
Yes: See the [Laboratory Glassware Waste Disposal](#) section.
No: Go to Step 6.

Step 6. Is the debris contaminated with a substance listed in the section [Liquid Non-Hazardous Chemical Waste Disposal](#)?
Yes: The contaminated debris can be disposed of in the regular trash.
No: Collect and store the contaminated debris as described in the waste container and storage guidelines listed below: dispose



Waste container for ethanol, acetone, methanol, isopropanol.




Waste container for mineral spirits.



Waste containers for chemicals used in NMR experiment

Follow Directly the Recommendations of Safety Working



Laboratory Sharps

Definition

Materials that qualify as “sharps” are defined at the state level and shall be disposed of as Potentially Infectious Medical Waste (PIMW). In Illinois, the Illinois Environmental Protection Agency (IEPA) has designated the following material (used or unused) as sharps:

- Any medical needles,
- Syringe barrels (with or without needle),
- Pasteur pipettes (glass),
- Scalpel and razor blades,
- Blood vials,
- Microscope slides and coverslips,
- Glassware contaminated with infectious agents.

NEVER dispose of these items in SDCs:

- Plastic items (except for syringes),
- Beverage containers (no pop cans!),
- Non-biologically contaminated laboratory glassware,
- Solvent/chemical bottles,
- Light bulbs,
- Any paper materials,
- Pipette tips,
- Plastic pipettes,
- Aerosol cans or cans of any type,
- Scintillation vials,
- Any item with liquid (except for blood in vacutainer tubes).



**Waste
container for
sharps**



Outline



V. Take a Lab tour ! It will be virtual tour.

VI. Let's get started
electronic logbooks
digital scopes



Outline

- I. Goals of the course
- II. Teamwork / grades / expectations from you
- III. Syllabus and schedule
- IV. Your working mode
In class and “after hours” access
Safety, Responsibility
Home and away computing
- V. Take a Lab tour !
- VI. Let's get started**
electronic logbooks
digital scopes



How to record the data

- Work together
- Write down the equipment used
- Make a diagram of the setup
- Note the settings of dials, switches, gauges
- Take a digital photo if appropriate (**we have prepared and will prepare more pictures of the setups equipment etc.)**
- Use a software drawing program to make a detailed sketch
(**PowerPoint works this very well**)



How to record the data

- Use the eLog (see next).
- Write down what you did in real sentences.
- Provide enough detail that you can reconstruct later what you did!
- How will you look at the data later?
- Do you have enough information?
- Did the equipment perform as expected?



How to record the data

- Many experiments require you to “change and measure” something by hand
 - Make a table in a paper logbook or put the data directly into electronic worksheet (*preferable*).
 - Make a “quick sketch” of your by plotting the data using OriginPro or other software

Looking on the graph you can answer the questions:


- Do you have enough points?
- Do you have any obvious anomalies?
- You can repeat points but do not throw them out.
Use other measurements to check reliability



How to record the data

- Many experiments have built-in, computer-based data acquisition (DAQ)

- You will not have time to fully understand the DAQ, but



- Be sure you know functionally what it is doing – ask
- A good idea is to make test measurements of something you know
- As before, anomalies? enough points? uncertainties?



Where to exchange, store and retrieve course information.

P403 Lab server

<\\engr-file-03\PHYINST\APL Courses\PHYCS403>



Connecting to the PHYS403 server

Connect to VPN following the instructions on the UIUC VPN website:

<https://techservices.illinois.edu/services/virtual-private-networking-vpn/download-and-set-up-the-vpn-client>

To connect to the PHYS403 Server:

- **Connect to the VPN first, then enter the following as the share to connect to:**
 - **Mac users:** Open Finder: Go: Connect to Server, type in address:
`smb://engr-file-03.engr.illinois.edu/PHYINST/APL Courses/PHYCS403`
 - **Windows users:** Open Windows Explorer, type in address:
[\\engr-file-03.engr.illinois.edu\PHYINST\APL Courses\PHYCS403](https://engr-file-03.engr.illinois.edu/PHYINST/APL Courses/PHYCS403)
- **When prompted for username and password, enter:**
“Uofl\[your netID]” and “[your netID password]”




Where to exchange, store and retrieve course information.

(i) Your data, projects, tables etc

\\engr-file-03\PHYINST\APL Courses\PHYCS403

There is a lot **useful** and not very useful stuff in many folders you can find there



**“Useful”
folders are
shown in red
frames**



Where to exchange, store and retrieve course information.



(i) Your data, projects, tables etc

\\engr-file-03\PHYINST\APL Courses\PHYCS403

Each student has a folder

- Library
- MuonData
- Students**
- Summer 2020 online
- Summer research


Store all experiment related materials in corresponding folder



Where to exchange, store and retrieve course information. (i)




Your data, projects, tables etc

An example of the “smart” structure of folders containing the raw data and data analysis projects



Where to retrieve course information.

Manuals, papers, setup diagrams and other useful materials




Manuals, papers, *setup diagrams* and other useful materials



α -range experiment setup diagram





Where to retrieve course information.

Setup diagrams – do not use cellphones to take the image of the setup from manual – for most setups we have PowerPoint projects with setups.



Where to retrieve course information.




Material Prepared for Online Teaching



“Journal club”

Lectures – Tuesday’s 3pm

Journal Club – Thursday’s 3pm



<http://ajp.aapt.org/#mainWithRight>

<http://www.nature.com/nature/index.htm>

<http://www.scientificamerican.com/>


<http://www.sciencemag.org/journals>



<http://publish.aps.org> or <http://prola.aps.org/>

“Journal club”

Walking with Coffee: Why Does it Spill?




Growth of Diamond Films from Tequila

J. Morales^{1,2}, L. M. Apátiga², V. M. Castaño²


1. Facultad de Ciencias Fisico Matemáticas, Universidad Autónoma de Nuevo León
2. Centro de física Aplicada y Tecnología Avanzada, Universidad Nacional Autónoma de México



Fabrication and Characterization of Ultrathin Three-Dimensional Thermal Cloak



(Credit: Guennea)



Student #1

University of Illinois at Urbana-Champaign

The Physics of Beer Tapping

PRESENTATION BY JOSEPH MIRABELLI

JAVIER RODRÍGUEZ-RODRÍGUEZ, 1,* ALMUDENA CASADO-CHACÓN, AND DANIEL FUSTER

1 FLUID MECHANICS GROUP, CARLOS III UNIVERSITY OF MADRID

2 CNRS, UNIVERSITÉ PIERRE ET MARIE CURIE

“Journal club”

Journal Access

If you cannot access journal papers using VPN, go to UIUC’s library proxy test site and enter the address of the paper you want to read:


<http://www.library.illinois.edu/proxy/test/>

Recommended journal websites

- **American Physical Society Journals:** <https://journals.aps.org/about>
- **Nature:** <http://www.nature.com/nature/index.html>
- **Science:** <http://www.sciencemag.org/journals>
- **American Journal of Physics:** <http://scitation.aip.org/content/aapt/journal/ajp>



Entering the e-Log ...



Home
Course Schedule
Gradebook
Course Description
Course Grading
Contact Information
Experiment Information
Lectures
Final Oral Session Abstracts
References
Online Materials
E-LOG
Section Information

PHYS 403 Summer 2020 **I**

Home page

Announcements


Link to e-Log

Welcome

Please see the [course description](#) for an explanation of how this course works. It may seem complicated at first, but all the pieces do work together to enhance understanding. Also, please consult the [schedule](#) to help you keep track of what is



Entering the e-Log ...



Please login

Username

Password

☒ Keep me logged in on this computer for the next 31 days or until I log out




[Register as new user](#)

Submit

Use your University Username and Password



Entering the e-Log ...

phys403-2020	
	Physics 403 Spring 2020  PHYS 403 Spring 2020 Semester
	Physics 403 Summer 2020  PHYS 403 Summer 2019 Semester
	Physics 403 Fall 2020  PHYS 403 Fall 2020 Semester



Entering the e-Log ...

Physics 403 Spring 2020 | Physics 403 Summer 2020 | Physics 403 Fall 2020

PHYS 403 Fall 2020 Semester

Navigation icons: [Back] [Forward] [New] [Reply] [Edit] [Find] [Delete] [Copy to] Phys403-Staff E-Logs

Message ID: 1 Entry time: 08/23/20 11:31

Author:	Eugene Colla
Experiment:	Intro
Post Type:	Other
Subject:	Welcome to Physics 403 course

Dear all,

Welcome to Modern Physics Lab course.

We will meet online on August 25th and next in person on August 27th

Eugene



e-logs: First a brief tour

How to use it

- **Pause and summarize your work at natural stopping points in the action. This is useful for particular findings and measurement sequences.**
- **Along the way, save data, plots, scope shots to your folder on the server.**
- **Near the end of the class, add a summary/conclusion, indicate future directions, and make sure the e-log provides a rather complete overview of the highlights of your work. Upload your plots, scope shots, etc. and describe the data.**

e-logs: Making a post ...

- **Create a New Post**
- **To create a new post, click "New" from the menu bar.**
- **Fill in the *Author, Experiment, Post Type, and Subject***

If the post is written by more than one person, use a comma separated list.

Be sure the Author name is the same you used when registering so that you can edit/delete the post if necessary.



e-logs: Making a post ...

Author:	Your name and your partner's name
Experiment:	General
Post Type:	How-To
Subject:	Day [#]: brief description of work

Goal: Be specific. Not, “Learn about experiment,” but, for example, “In helium below temperatures of 2.17K, a second sound due to thermal effects becomes measurable. We will measure second sound using a resonant cavity...”

Settings / Equipment Notes: Note important environmental and experimental parameters such as atmospheric pressure, settings on equipment, etc.

[Time Range 1]: Give time range, not just “before tea.”


- Note important steps and results
- Include plots, photos, or scope shots in attachments below
- Use bullet points to make it easy to read

[Time Range 2]: ...

Conclusions & Future Plans: What did you find and what is the next step? Be specific. Not, “We measured decay times,” but, for example, “Ruby #2 sample with higher concentration chromium was observed to decay with a form...”

e-logs: Making a post ...

Elog records should contain the information about parameters of the experiment and that is why we suggest you to use the templates (<\\engr-file-03\PHYINST\APL Courses\PHYCS403\Common\elog templates>):



Ferro1	1/15/2020 2:39 PM	Microsoft Word D...	17 KB
Ferro2	1/15/2020 2:05 PM	Microsoft Word D...	15 KB
Ferro3	1/15/2020 2:45 PM	Microsoft Word D...	15 KB
Superconductivity	1/19/2020 8:58 PM	Microsoft Word D...	16 KB
Superconductivity_mutual inductance	1/19/2020 9:01 PM	Microsoft Word D...	15 KB
Tunneling	1/19/2020 9:07 PM	Microsoft Word D...	15 KB



e-logs: Making a post ...

Copy and Paste the template (table) into the record and fill it up with numbers corresponding experiment parameters

Message ID: 365 Entry time: 01/14/20 16:34

Author:	Eugene Colla
Experiment:	Ferroelectric (Dielectric)
Post Type:	Measurement
Subject:	example of using of the template

BaTiO ₃		BT1		Sample area: 4.01 mm ²		Sample thickness: 0.8 mm
File name	Folder	T range (K)	Frequency (Hz)	V _{AC} (V)	V _{DC} (V)	Comments
14JAN20_s1	Data:student:BTO:set1	300-100K				



Some General Physics 403 Rules.



No Food or Drinks in Lab

