


Electrical Cross-Correlation Spectroscopy: Measuring Picoliter-per-Minute Flows in Nanochannels

Klaus Mathwig, Dileep Mampallil, Shuo Kang, and Serge G. Lemay, Phys. Rev. Lett. 109, 11 (2012).



Angela Chen

Li Chen

Patrick Coleman

Karl Decker

University of Illinois

November 30, 2012

Outline

- Context for experiment
- Experiment and analysis
- Moving beyond the experiment



Outline

- Context for experiment
- Experiment and analysis
- Moving beyond the experiment



Prior Work

- Parallel flow control (PFC) system for syringe driven nanofluidics—measured nanochannel flow rate at the value of 30 pL/min.
- Fluorescence correlation spectroscopy—by measuring the rates of decay of spontaneous concentration fluctuation, we can determine chemical kinetic constants and diffusion coefficients.

Why use cross-correlation spectroscopy?

- Estimate ultralow liquid flow rates of electrochemically active molecules.
- Realize multiple detectors in more complex nanochannel networks.
- Avoid using bulky optical equipment (lab on a chip compatibility).

Outline



- Context for experiment
- **Experiment and analysis**
- Moving beyond the experiment



Time-of-flight measurements

- Detect current in nanogap transducers.
- Perform time-of-flight measurement of fluctuations using cross-correlation functions.


Parallel Flow Device



Schematic of parallel flow. Mathwig et al., 2012


- Slow flow rates: 5 $\mu\text{L}/\text{h}$ to 50 $\mu\text{L}/\text{h}$
- Microchannel bonded to nanochannel.

Nanogap transducers measure fluid flow





- 130 nm channel between electrodes.
- Repeated oxidation and reduction of molecules.
 - Enhanced current
 - Current proportional to number of molecules

Two detection volumes track fluctuation




- Two detection electrodes.
- Fluctuation measured by both electrodes as it flows downstream.

Cross-correlated current signals yield flow rates




Peaks indicate time of flight (τ)




- Higher flow rates correspond with lower time of flight.
- Noise at $\tau > 1$ s is due to finite length of signals.

Simulation agrees with experiment




Experiment




Simulation

Auto-correlation yields time-of-flight




Experiment



Simulation

Nanochannel flow rate varies linearly with pump flow rate



- Adjusted data accounts for adsorption.
- Remaining disparity blamed on bulging microchannel.

Outline



- Context for experiment
- Experiment and analysis
- **Moving beyond the experiment**



Huge improvements in tiny flow rate detection

- Flow rate measurements of 10 pL/min – slowest rate ever reliably measured
- Limited by signal-to-noise issues
- Possible improvements: “Improved instrumentation or reduced channel height”


Innovating towards...




<http://web2.clarkson.edu/projects/nanobird/2.4.php>

- Compatible with complete “Lab-on-a-chip”
- Ability to identify single-molecules


Summary



Nanogap transducers



Small flow rate detection



Lab-on-a-chip applications

- Since the article was published on September 12, 2012, no papers have cited it yet.

Questions?

